



Fast APM Inflation-Deflation Pump Kit



Inflate and deflate alternating pressure mattresses (APM) quicker, reducing setup and pack down time.

- Inflate or deflate alternating pressure mattress in under 2 minutes¹.
- Portable and light weight.
- Mains powered pump with high volume output for rapid results.
- Includes adaptors to suit all current Novis mattresses².
- Adaptors made from durable reinforced silicone rubber suitable for heavy duty usage.
- Especially useful for maintenance, ward and nursing staff.



1 Will fully inflate all cells of the Premium 5/8 APMs and equivalent. The three static head cells of the Premium 5/8 and equivalent will require manual deflation by opening CPR valve.

2 Except Tridien Sentech APMs/turning mattresses.

CODE	APMFIDK
INFLATION/DEFLATION PUMP	1x 240v mains powered high volume pump
ADAPTORS SUITABLE FOR (NOVIS MODELS)	ProCair ProCair Plus ProCair Prime ProCair Pro BetterLiving APMs Premium 10 Premium 9 Premium 8 Premium 5
WARRANTY	1 year